



Endurus Silica's exceptional consistency is the foundation of Endurus column performance. Stringent quality control standards produce unmatched batch-to-batch reproducibility. Column packing technology minimized column-to-column variation ensure high efficiency, excellent peak symmetry and extremely low operating pressure.

Capacity Factors for Endurus HP C18 Columns

◆ Biphenyl	■ Toluene	▲ N-butylbenzoic acid
● Amitriptyline	× Nitrobenzene	

Consistent batch-to-batch capacity for acidic, basic & neutral compounds



Test Conditions: Endurus C18 (250 x 4.6 mm, 5 μ m), Eluent: CH₃CN - 0.05M KH₂PO₄, pH 3.2 (60:40), Flow rate: 1 ml/min, Detection: 254 nm